


CARSON

DR-200 ezRead Electronic Reading Aid


1. Installation of Battery



2. Using Power Adapter (Alternative Power Source)




3. Connecting the Digital Reader to Television Set



200mm

4. Using the Digital Reader



Specifications:

- Magnification: 8X ON 21" TV
- Field of View: 56X42mm
- Illumination: Built-in LED
- Resolution: 250K PIXELS
- Interface: TV NTSC
- Power: 3X AAA Batteries(Not included)
110V Power Adapter

Warning

To avoid possible permanent eye damage, do not place lighted Digital Reader directly upon human eye.

MADE IN CHINA

2037190001